



Get a good start
with your Kirkland
Signature 6.0
Hearing Instruments
with the
Choice App



My Information

Please provide the following information and bring this brochure with you to your next appointment:

My smartphone is:

My Apple / Google Play ID is: _____

Password: _____

My other streaming devices:

Welcome to a new generation of hearing aids

Your Costco Hearing Aid Center professional has placed your order.

You will soon be ready to experience our Kirkland Signature 6.0 hearing instruments, including the ability to stream sound directly from your device to your hearing instruments*. In order to prepare for your next appointment, you are kindly asked to download the Choice app.



*Direct streaming available with Apple devices. Phone Clip+ required for direct streaming from compatible Android devices

How to download the Choice app:

Step 1: Make sure your smartphone is compatible with your Kirkland Signature 6.0 hearing aids.

For a current list of compatible devices, consult <http://ks6userportal.com/compatibility/>

Step 2: Download the Choice app.

From the App Store (Apple devices)

1. Tap the App Store icon. Select **iPhone only** setting. Please note that this will allow you to download to an iPad and an iPod touch too. This setting is only to locate the Choice app.
2. Type **Kirkland Signature Choice** in the search field.
3. Once the Choice app screen is displayed, tap **Free** (It will then change to **Install**).
4. Tap **Install** to load the app on your device.
Note: Although this is a free app, you will be prompted to sign into your iTunes® account.
5. Type your Apple ID and password to complete the download process.

Step 3:

From Google Play (select Android devices)

1. Go to Play Store on your **Android phone** and search for Kirkland Signature Choice
2. Tap on **Install**, then **Accept** so the app can access certain features on your phone such as Bluetooth
3. Your phone will now install the app
4. When installation is complete, tap on **Open** to start up your Kirkland Signature Choice app

You are now ready to pair your hearing instruments with your smartphone.

If you need assistance with your smartphone, please contact your local Apple store or smartphone provider to make an appointment.





Distributed by:

Costco Wholesale Corporation

P.O. Box 34535
Seattle, WA 98124-1535
USA
1-800-774-2678

Costco Wholesale Canada Ltd

415 W. Hunt Club Road
Ottawa, Ontario
K2E 1C5 Canada
1-800-463-3783

Find out more at www.costco.com

www.costco.ca

Kirkland Signature 6.0 Hearing Instruments are compatible with iPhone 6, iPhone 6 Plus, iPhone 5s, iPhone 5c, iPhone 5, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini 2, iPad mini and iPod touch (5th generation) using iOS 7.X or later. Apple, the Apple logo, iPhone, iPad, iPod touch and iTunes are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License.



MK604584

MK604584 Rev B 2015.10