



Making the most
of your Kirkland
Signature 6.0
Hearing Instruments
with the
Navigator app



My Information

Please provide the following information and bring this brochure with you to your next appointment:

My smartphone is:

My Apple / Google Play ID is:_____

Password: _____

My other streaming devices:

Welcome to a new generation of hearing aids

Your Costco Hearing Aid Center professional has placed your order. You will soon be ready to experience our Android-compatible hearing instruments*, including the ability to stream sound from your Android device to your hearing instruments. In order to prepare for your next appointment, you are kindly asked to download the Navigator app.



*The optional Phone Clip+ accessory is required to use the app.

How to download the Navigator app:

The Navigator app is free to download.

On your mobile device:

- Go to Google Play and search for **Kirkland Signature Navigator** app.
- Click “**install**” and download directly to your mobile device.
- “**Accept**” the request to access your Bluetooth connection information.

If you do not have an account on Google Play, please follow the Android instructions on how to set up an account when entering Google Play on your mobile device.

The Navigator app is compatible with most Bluetooth-enabled devices.


How to pair your Kirkland Signature 6.0 hearing instruments to the Phone Clip+ and an Android phone

Step 1: Pair mobile device to the Phone Clip+:

1. Make sure the Phone Clip+ is charged. Turn it on and remove the silver cap.
2. Place your mobile device next to the Phone Clip+ and make sure the mobile device is on.
3. Find the connectivity settings on your mobile device and make sure Bluetooth is on.
4. With the tip of a pen or similar object, press the small blue pairing button on the back of the Phone Clip+ one time. This will put the Phone Clip+ in pairing mode.
5. On the mobile device, choose to search for Bluetooth devices.
6. Your mobile device will display a list of Bluetooth devices discovered. Choose **Hearing Aid Phone** from this list.
7. If you are asked to enter a passcode, enter “**0000**” (four zeros).
8. Your Phone Clip+ should now be paired to your mobile device. Some mobile phones may ask which Bluetooth service you wish to enable. Select **Headset**.

Step 2: Pair the Phone Clip+ to your Kirkland Signature 6.0 hearing instruments:

1. Remove the silver cap from the Phone Clip+.
2. With the tip of a pen or similar object, press the white pairing button on the Phone Clip+.
3. Within 20 seconds, open and close the battery doors of your hearing instruments.



How to start using the Navigator app

After downloading the app and pairing it to both the hearing instrument and Phone Clip+, simply **tap** the launch icon on your home screen to activate the app.

1. You will be asked: Do you want to connect now or run in demo mode? Choose **“connect now”**.
2. Terms of use: You will be asked to familiarize yourself with the printed user guide for the Phone Clip+. It is necessary to **accept** the terms of use to proceed.
3. We recommend that you disable any lock sounds and keyboard clicks on your mobile device. This is usually done under **“Settings”** and **“Sound”**. Lock sounds and keyboard clicks may interfere with Phone Clip+/app operation.
4. You will be directed to find tutorials under **“Settings”** to demonstrate the functionality of the app and help you pair your mobile device with your Phone Clip+.



Distributed by:

Costco Wholesale Corporation

P.O. Box 34535
Seattle, WA 98124-1535
USA
1-800-774-2678

Costco Wholesale Canada Ltd

415 W. Hunt Club Road
Ottawa, Ontario
K2E 1C5 Canada
1-800-463-3783

Find out more at www.costco.com

www.costco.ca

Android and Google Play are registered trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.



MK604585 Rev B 2015.10